

# Brillant 200 Cut Off Machine

10  
inch



## Simple. Spacious.

The **Brillant 200 Cut Off Machine** is the ideal selection for a basic, wide capacity metallurgical saw.

Constructed with a solid welded aluminum frame and stainless steel lift door, the Brillant 200 is tough and corrosion resistant. A 10 inch cut off wheel and a 5 HP motor provide a cutting capacity of 3.5 inch diameter cross section.

The large retractable door provides optimum access

and the optional illuminated interior provides an excellent view of the cutting area.

The Brillant 200's Vertical Clamping System provides quick and versatile left and right handed options. Openings on both sides of the cutting chamber allow for cutting of long, continuous parts.

The Brillant 200 offers labs with a smaller sample volume a superb choice of a manual cut off machine.



# Brillant 200 Cut Off Machine

Max Cutting Capacity:  
90 mm / 3.5" Diameter

Max Cut Off Wheel Size:  
256 mm / 10"

Arbor Size:  
32 mm / 1 1/4"

Table Size:  
350 x 288 mm / 13.8 x 11.4"  
12 mm T-slot

Motor: 5 HP

Power:  
220, 480 V / 60 Hz,  
3 Phase

Dimensions W x H x D:  
32 x 29 x 26"

Weight: 166 lbs

### Optional:

Recirculation System:  
45 L / 12 Gallons

SystemLab Base Cabinet  
One door, 5 drawers  
Dimensions W x H x D:  
35 x 32 x 32"



Wide open work area incorporates a **large T-slot Table** with replaceable stainless steel top plates



The optional **LED task light** is easy to position to provide optimal illumination



Unique, **Dual Chamber Recirculation System** with stainless steel basket and filter bag preserves coolant longer



The **BR200 Vertical Clamping System** is integrated into the stainless steel table for a quick-clamping option

