

Topas 130 Dosing System

Automatic



Simple. Integrated.

The **Topas 130 Automatic Dosing System** makes dosing diamond abrasive easy. With space for three diamond suspensions and one lubricant, each solution can be adjusted to dose the perfect amount

for proper preparation of your samples.

While the **Topas 130** is designed to integrate flawlessly with ATM's line of grinder/polishers, it can also be used with any other brand.

Topas 130 Automatic Dosing System

Variable Dosing:
Manual or Automatic
using Peristaltic Pumps

Dosing Intervals:
2 seconds - 5 Minutes

Number of nozzles:
3 Suspensions
1 Lubricant

Power:
110 V / 60 Hz, 1 Phase

Dimensions W x H x D:
7 x 10 x 8"

Weight: 5 lbs



Individually control the
precise, variable dosing of
suspensions and lubricants
with two second to five
minute dosing intervals



Exact droplet dosing with
reverse rinsing function
keeps lines clear and ready
for use



Mager's line of high
performance diamond
suspensions are made
specifically for manual and
automatic dosing systems

